# LTS2000 with ISM50 Interactive Sensing Module

© Harrith M. Hasson, MD



#### Introduction

•The LTS2000 is a physical training simulator for developing and testing laparoscopic surgical skills.

•The ISM50 Interactive Sensing Module adds electronic scoring and feedback to the existing LTS2000 platform.

## **LTS2000 Physical Simulator**

 Enclosure with multilayered abdominal wall

 Adjustable floormat to which various exercise modules are attached

 Camcorder or video laparoscope



## LTS2000 with ISM50 Upgrade

- Carousel with physical sensors
- System controller that links carousel to camera, monitor and PC; overlays score data on video monitor
- Benefits:
  - Real-time user feedback
  - Quantification of performance metrics
  - Facilitates data storage and analysis



## Sensor Carousel

- Incorporates 5 laparoscopy training tasks on one platform
- Potentiometer switches tasks as carousel is rotated



#### System Controller

 Controls task timing and error logging and sends data to PC



 Displays progress and score on video monitor (overlay)



#### **Practice Mode**

 "Personal Best" achievement

 Data displayed but not recorded

PC not required




## **Testing Mode**

 PC software shows real-time feedback

- Data logged to Excel file:
  - User ID / experience level
  - Test data including results
  - Cumulative scores
  - User login/logout





## **Coordination Tasks**


**Peg insertion** Insert 6 pegs into 9 holes






**Ring manipulation** Guide 2 o-rings down wire guides, avoiding contact





Cannulation Thread flexible rod through narrow tube and remove






# Knot Tying Tasks


**Knot tying** Cinch knots around vertical and horizontal tubes





Knot integrity Tie knot between two plates; test with 2.5 lbs tensile force





#### Conclusion

- Virtual Reality systems provide excellent visual representation but lack realistic physical sensations due to limited haptic feedback
- LTS2000 with ISM50 is the first simulator to combine:
  - Realistic physical experience
  - Real-time interaction
  - Computerized data storage

 LTS2000 training complements both laboratory and VR simulator training

### **Further Information**

Harrith M. Hasson hmh@lts2000.com www.lts2000.com